

# INSTRUCTION MANUAL FOR PLANK'S CONSTANT APPARATUS

Plank's constant apparatus has been designed to determine the value of Plank's constant 'h' by a photo cell.

## Apparatus Used :

Photo emissive cell mounted in a box provided with a wide slit D.C. Power Supply, set of filters, filter stand, light source etc.

## Formula Used :

$$h = \frac{e (V_2 - V_1) \lambda_1 \lambda_2}{c(\lambda_1 - \lambda_2)}$$




Fig - 1




Fig - 2

e	=	Electronic charge
$V_2$	=	Stopping potential * corresponding to wave length 2
$V_1$	=	Stopping potential * corresponding to wave length 1
c	=	Velocity of light.

\* Minimum negative potential applied to anode to reduce the photo electric current to zero.

## For Plank's Constant :

- (i) Keep the left hand side switch on the panel towards sensitive side and right hand side switch towards 1.0 V side. Switch on the unit. Now set the  $\mu A$  reading to zero with the help of potentiometer marked with zero adjust.
- (ii) The circuit connections are made as shown in diagram (Fig - 1). Be careful about the polarity shown in diagram.
- (iii) A light source is arranged. The light is allowed to fall on the tube. The distance between tube and light source is adjusted such that there is a deflection of about 8 to 10 div. in  $\mu A$ . Now a suitable filter (say green) of known wave length is placed in the path of light (in the slit provided) say it is with wave length  $\lambda_2$ .



- (iv) A deflection is observed in the micro-ammeter. This deflection corresponds to the zero anode potential.
- (v) A small -ve potential is applied on the anode. This voltage is recorded with the help of voltmeter provided (1.0 volts range).
- (vi) The negative anode potential is gradually increased in steps and each time corresponding deflection is noted till the micro-ammeter deflection reduces to zero and this is stopping potential  $V_2$  corresponding to filter with wave length  $\lambda_2$ .
- (vii) The experiment is repeated after replacing the green filter with blue and red filters. Say with wave length  $\lambda_2$  and  $\lambda_3$  respectively and stopping potential  $V_1$  and  $V_3$  are noted.
- (viii) Taking negative anode potential on x-axis and corresponding deflections in micro-ammeter on y-axis, graphs are plotted for different filters.
- (ix) By using above values Plank's Constant 'h' is calculated by the formula given. Standard values of e, c and wave length of standard filters are given below.

$$e = 1.6 \times 10^{-19} \text{ coulombs}$$

$$c = 3 \times 10^8 \text{ m/sec.} = 3 \times 10^{10} \text{ cm/sec}$$

$$\text{Wave length of green filter } \lambda_3 = 5645 \times 10^{-10} \pm 2\% \text{ meter.}$$

$$\text{Wave length of red filter } \lambda_1 = 6143 \times 10^{-10} \pm 2\% \text{ meter.}$$

$$\text{Wave length of blue filter } \lambda_2 = 5265 \times 10^{-10} \pm 2\% \text{ meter.}$$

### ***For Photo Cell Characteristics :***

- (i) Keep the left hand side switch towards normal and right hand side switch towards 50 volts side.
- (ii) The circuit connections are made as shown in the diagram (fig - 2).
- (iii) Light is allowed to fall on the tube.
- (iv) A small +ve potential is applied on the anode and corresponding reading of micro ammeter is noted.
- (v) The potential is gradually increased in steps and each time corresponding reading in Micro ammeter is noted.
- (vi) A graph is plotted by taking V on x-axis and  $\mu\text{A}$  on y-axis.